English

Italian/Italiano

German/Deutsch

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Swedish/Svenska

EN

## INSTRUCTIONS FOR USE

### MODELLATION

Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

### SPRUIING

Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or in-direct technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

### INVESTMENT

Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

### BURN-OUT

The suggested burnout temperature: 800-850°C/1470-1560°F

### MELTING AND CASTING

Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Always use fresh alloy. Depending on the type of casting machine, follow the manufacturers instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature. Torch and induction casting: The ingots slump down, look for oxide break and cast within next 2-4 seconds.

Casting Temperature: 1450-1500°C/2640-2730°F

### METAL PREPARATION

Carefully divest and clean the object with Al<sub>2</sub>O<sub>3</sub>. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Blast surface with 50-100 micron Al<sub>2</sub>O<sub>3</sub> at max. 4.5 bar/65 psi pressure. Subsequently, clean with ultrasonic or steam. Avoid inhalation of dust during grinding!

### OXIDATION

Not required.

Quality check: recommended to check the alloy surface.

Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650°C/1200°F and increase the furnace temperature to 950°C/1740°F under full vacuum with 1 min. hold time at top temperature. Use ceramics according to ISO 9693 with firing temperature of maximum 1000°C.

### SOLDERS AND FLUXES

Design the soldering patty as small as possible and preheat it in the furnace at approximately 600°C/1112°F. The soldering patty should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Pre-Solder: Special High Fusing White Ceramic Solder, High Fusing Bondal Flux NP

Post-Solder: Low Fusing White Gold Solder, High Fusing Bondal Flux NP

Laser weld material: Laser Ceramic White

### INDICATIONS

Recommended for crowns, telescope crowns, conus crowns, posts, short- and long span bridges, PFM crowns, implant superstructures, partial dentures.

CONTRAINDICATIONS  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

### SIDE EFFECTS

In individual cases, sensitivity or allergies to components of this alloy may occur.

### INTERACTIONS

Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

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